


STIM onTrack™
Mobile App



Connecting the STIM onTrack™ App to your Bone Growth Therapy Device



In order for the STIM onTrack App to connect to your bone growth therapy device, the device must be on and providing a treatment.

1. The first time you launch (open) the STIM onTrack App, you will be prompted to 'Connect to Your Device'.
2. Scan the serial number barcode on your bone growth therapy device by pressing the button 'Scan Device Barcode' or by manually typing in the numbers in the open text box.
3. Push the 'Connect to Device' button. The appearance of the Treatment Screen is an indication that the device has paired with the app successfully.

⚠ Note: The Orthofix STIM onTrack Mobile App is supported on iOS devices. The bone growth therapy device will not communicate with any other device or app other than the official Orthofix STIM onTrack App from Orthofix available on the App Store.

⚠ Note: Bluetooth on your Apple device must be turned on.